



# ELAN

Straight Indoor Stair Lift

 **BRUNO****MADE IN USA**



# BRUNO ELAN

#1 Selling Stairlift Made in the USA

## Easy to Use

The Elan smoothly glides up and down your staircase with the press of a single button.


## Space Saving

The Elan's purpose-built vertical rail maximizes open space on steps for family, friends and pets. The armrests, seat and footrest fold up, creating a compact storage profile.

## Dependable

The Elan is battery-powered and continuously charges anywhere it stops along the rail, delivering reliable operation even in the event of a power outage.






## AVAILABLE OPTIONS

Make rotating the Elan seat at the top of the steps automatic with the selection of the **power swivel** option.

Easily fold the footrest up into the storage position by choosing the popular **power footrest** option.

Maintain open space at the bottom of the staircase with the addition of a manual or **power folding rail** option.

*Not all options shown.*





## TECH SPECS

**Lift Capacity:**

300 lb (136 kg)

**Seat Size:**

14.1 in (358 mm) deep x  
15.38 in (391 mm) wide

**Width Between Armrests:**

20.5 in (520 mm)

**Footrest Size:**

10.5 in long x 10.75 in wide  
(267 mm x 276 mm)

**Vertical Rail:**

Installs within 5 in (127 mm) of the wall

## WARRANTY

**Limited Lifetime Warranty on Major**

**Components** (motor, gearbox, drivetrain and stationary rail system) for the lifetime of the Bruno stair lift as long as the original purchaser owns the stair lift.

**2-Year Warranty** on all parts after date of installation by an Authorized Bruno Dealer.

**30-Days Labor Coverage** after date of installation by an Authorized Bruno Dealer for any manufacturing defect.



Scan for more  
info & videos



Bruno's mission is to help make life easier for people with mobility challenges by delivering the highest quality stair lifts, scooter lifts, platform lifts and home elevators.

© Bruno Independent Living Aids, Inc. 2026  
2026-05-05-elan-mkt-sl-6.indd